


APRIL - JUNE 2019 DURALAST 'LITE'




MADE IN UK

The Next Generation of Hand Instruments

- ✓ Good instruments do not need to be expensive.
- ✓ Duralast 'Lite' means less hand fatigue.
- ✓ Laser Etched for ease of identification.
- ✓ High carbide steel produces razor sharp long lasting edges.

**BUY 4
GET 1
FREE!**



Heat treated tips
for maximum
hardness & razor
sharp edges

Duralast 'Lite'
weighs just 13 gm

Peaks & trough
profile provides
secure low
pressure grip

\$34.90

Each



Unit 4, 130 - 132 Bayfield Road East
Bayswater North
VIC 3153
P: 1800 220 211 | F: 9738 0822
sales@matrixdental.com.au


New Light Weight Handle

- Feather light 13gm instrument
- Large 10mm diameter handle.
- Machined to give tactile control & less finger fatigue.
- Easy to clean surfaces.


New Cutting Edge Metallurgy

- Manufactured from 440C grade stainless steel.
- Heat treated to high Rockwell hardness.
- Polished & sharpened to a razor edge.


CURETTES




27200N
Curette Gracey 1-2




27201N
Curette Gracey 3-4



27202N
Curette Gracey 5-6




27203N
Curette Gracey 7-8




27204N
Curette Gracey 9-10




27206N
Curette Gracey 11-12




27208N
Curette Gracey 13-14



27250N
Curette Columbia 2L-2R



27251N
Curette Columbia 4L-4R



27252N
Curette Columbia 13-14

**BUY 4
GET 1
FREE!**


\$34.90
Each

ENDODONTICS

ENDO NEEDLES



ENDO NEEDLE
OPEN END



ENDO NEEDLE
CLOSED END



25027
Endo Needle
Closed End
Side Vent 27G
\$9.00
Bag - 20

25030
Endo Needle
Closed End
Side Vent 30G
\$9.00
Bag - 20

25127
Endo Needle
Open End Side
Vent 27G
\$9.00
Bag - 20

25130
Endo Needle
Open End
Side Vent 30G
\$9.00
Bag - 20


IRRIGATING SYRINGES



49001
MGuard Irrigating Syringe

\$33.00
Pkt - 50

DENTAL DAMS NON-LATEX



68100
MGuard Non-Latex Dental Dam
15cm x 15cm Blue

\$29.00
Pkt - 20

- Holes can be placed close together as the material readily stretches.
- It is extremely flexible and doesn't tear when stretching it.
- Easy to place and floss it through without the need for lubrication of the dam.

LEDERMIX



007698
Ledermix Paste

\$139.00
1x5g



007699
Ledermix Cement Powder

\$99.00
1x3g



007696
Ledermix Hardener

\$28.00
1x2.5mL



007693
Dropper + Mixing Pad

\$3.00
Each

K FILES

Blade grinded with a square section then twisted in order to obtain cutting power and flexibility. Helix angle higher than on reamer. Packs of 6.



Endo Files K - 21mm

Ø ISO	6	8	10	15	20	25	30	35	40	45	50	55	60	70	80
Colour															
Code	80306	80308	80310	80315	80320	80325	80330	80335	80340	80345	80350	80355	80360	80370	80380

Endo Files K - 25mm

Ø ISO	6	8	10	15	20	25	30	35	40	45	50	55	60	70	80
Colour															
Code	80406	80408	80410	80415	80420	80425	80430	80435	80440	80445	80450	80455	80460	80470	80480

H FILES

This instrument is directly grinded. High cutting efficiency, use in filing action. Packs of 6.



\$9.90
Pkt - 6

Endo Files Hedstroem - 21 mm

Ø ISO	8	10	15	20	25	30	35	40	45	50	55	60	70	80
Colour														
Code	80008	80010	80015	80020	80025	80030	80035	80040	80045	80050	80055	80060	80070	80080


Endo Files Hedstroem - 25mm

Ø ISO	8	10	15	20	25	30	35	40	45	50	55	60	70	80
Colour														
Code	80108	80110	80115	80120	80125	80130	80135	80140	80145	80150	80155	80160	80170	80180

\$9.90
Pkt - 6

PROTECTIVE APPAREL

SURGICAL MASKS



MGuard Level 2 Antifog Mask

- Meets Australian Standards AS 4381-2015
- 3 layer masks for appropriate fluid protection
- Large non-collapsing breathing chamber for greater comfort
- Correct overall sizing provides more secure facial fit



MGuard Level 3 Fog Free Mask "Protects You from aerolised microbes"

- Meets Australian Standards AS 4381-2015
- 4 layer masks for appropriate fluid protection
- Large non-collapsing breathing chamber for greater comfort
- Soft smooth cellulose inner lining
- Correct overall sizing provides more secure facial fit
- More protection for longer procedures
- Most masks are rated for only 20 minutes.

"High Fluid Risks Need High Level Protection"

**BUY 5
GET 1
FREE!**
LEVEL 2

70208
MGuard Mask 120
Barrier Level 2
Earloop Antifog

\$9.90

Box - 50

70206
MGuard Mask 120
Barrier Level 2
Tie On Antifog

\$9.90

Box - 50

70214
MGuard Mask 160
Barrier Level 3
Earloop

\$14.90

Box - 50

UVEX

EYEWEAR ANTIFOG SCRATCH RESISTANT



9192302

UVEX Pheos Hardcoat Scratch Resistant Antifog Lens Clear Adult



9192304

UVEX Pheos Hardcoat Scratch Resistant Antifog Lens Tinted Adult

PHEOS

- Antifog hard-coated lens
- Side protection
- Soft nose bridge

\$23.90

Each

**BUY 3
GET 1
FREE!**

EYEWEAR ADJUSTABLE SIDE ARMS



9195275

UVEX Skyper Black Frame Clear Supravision Lens



9195075

UVEX Skyper Black Frame H/C Clear Lens



9196065

UVEX Skyper Child Black Frame H/C Clear Lens

SKYPER

- Integrated side protection
- 5 point side arm adjustment
- Pressure free comfort fit

\$29.00

Each

**BUY 3
GET 1
FREE!**

EYEWEAR ANTIFOG OVER THE TOP



9169945

UVEX Super OTG Frame Clear Lens



9169946

UVEX Super OTG Frame Tinted Lens

OTG

- Fits over prescription spectacles
- Integrated brow guard
- 4 point length adjustment

\$29.00

Each

MGuard

**BUY 5
GET 1
FREE!**



25201

Prophy Cup Snap Blue Non Latex



25100

Prophy Cup Snap White Non Latex

\$13.90

Pkt - 100

**BUY 5
GET 1
FREE!**



25304

Prophy Cup Screw In Blue Non Latex



26410

Prophy Cup Screw In White Non Latex

\$13.90

Pkt - 100

**BUY 5
GET 1
FREE!**



26412

Prophy Cup Latch Blue RA Non Latex



26413

Prophy Cup Latch White RA Non Latex

\$13.90

Pkt - 100

**BUY 5
GET 1
FREE!**



75100

Prophy Brush Jcup RA White

\$39.90

Pkt - 100

WHILE STOCKS LAST, ALL PRICES INCLUDE GST, WE RESERVE THE RIGHT TO AMEND ERRORS & PRICES.

DENTAGENIE

PREVENTIVE

“87 % Floss Infrequently or Not at all”

HYGIENISTS TELL US PATIENTS HATE FLOSSING AND HAS A LOW RATE OF COMPLIANCE.

According to the American Dental Association (ADA), “the best way to remove decay-causing plaque is by brushing and cleaning between your teeth every day. Brushing removes plaque from tooth surfaces... cleaning between your teeth with floss or interdental cleaners removes plaque between your teeth where toothbrushes can’t reach.”*

Encourage your non-flossing patients to clean between their teeth by giving a low cost alternative.



LOW COST ENTRY

Small cylindrical shape fits spaces as small as 1.0mm

Size	Code	Description
SIZE 0	10060	For very fine tooth gaps
SIZE 1	10061	For small tooth gaps
SIZE 2	10062	For small to medium tooth gaps
SIZE 3	10063	For medium tooth gaps
SIZE 4	10064	For medium to large tooth gaps
SIZE 5	10065	For large tooth gaps
	10066	Professional Sampler Pack

*ADA News Release, American Dental Association, September 2000. Source: 199 Gallupp Study

\$0.20
Per Piece



\$29.00
Box - 12 Pkts

FLOSS AND TAPE



33021
Reach Dentotape 20m

\$22.90
Box - 6



33028
Reach Waxed Floss 50m

\$16.90
Box - 6



33040
Reach Flosser Cleanpaste Starter

\$34.90
Box - 6



33041
Reach Flosser Cleanpaste Refill

\$44.90
Box - 6

WHILE STOCKS LAST, ALL PRICES INCLUDE GST, WE RESERVE THE RIGHT TO AMEND ERRORS & PRICES.



UNBEATABLE PRICES

AS LOW AS \$11.00 PER BAG

ALGINATE



90100 Hydrogum 5 Alginat

\$12.00
1x453g



90116 Orthoprint Alginat

\$12.00
1x500g



90106 Phase Plus Alginat

\$11.00
1x453g



90118 Neocolloid Alginat

\$14.00
1x500g

Hydrogum 5 Alginat

- 5 days of dimensional stability
- Higher consistency
- Smooth, even impressions

Orthoprint Alginat

- Highly elastic, even with braces
- Pleasant flavour, reducing the emetic reflex
- Excellent impression readability

Phase Plus Alginat

- Colour change alginat
- Creamy consistency
- High final hardness

Neocolloid Alginat

- High Precision
- Extended working time 2 minutes
- Ideal for chromes

ALGINATE REPLACEMENT



Freealgin Alginat Replacement

- Cost effective alternative to alginates
- Higher accuracy than alginates
- Can be poured multiple times - up to 14 days stability
- Plus four tips

90151 Freealgin Alginat Replacement

\$49.00
Box - 2x50 mL

ELITE HD+ VPS MATERIAL



90300 Elite HD+ Putty Normal Set

\$120.00
Each



90301 Elite HD+ Putty Fast Set

\$120.00
Each



90200 Elite HD+ Tray Material Fast Set

\$55.00
Each



90206 Elite HD+ Light Body Fast Set

\$55.00
Each



90207 Elite HD+ Light Body Normal Set

\$55.00
Each



90204 Elite HD+ Monophase

\$55.00
Each



90205 Elite HD+ Medium Body Normal Set

\$55.00
Each

WHILE STOCKS LAST, ALL PRICES INCLUDE GST, WE RESERVE THE RIGHT TO AMEND ERRORS & PRICES.